


R6880D

R6870D

155 mm, 6"




R6890D




WILLIGER  
R6900D

165 mm, 6 1/2"



WILLIGER  
R6910D

170 mm, 6 3/4"



FREER-YASARGIL  
R6920D

185 mm, 7 1/4"




HALLE  
R6930D

180 mm, 7"




SYME  
R6940D

185 mm, 7 1/4"




FREER  
R6950D


185 mm, 7 1/4"




**R6960D**  
**UM309K**  
155 mm, 6"




**R6970D**  
**DU510K**  
140 mm, 5 1/2"



**R6980D**  
**FK300K - FK303K**  
160 mm, 6 1/4"




**R6990D**  
**FK301K - FK302K**  
165 mm, 6 1/2"




**R7000D**  
**OBWEGESER**  
**DU530K**  
170 mm, 6 3/4"



**R7010D**  
**DU531K - DU532K**  
200 mm, 8"



**R7030D**  
**DU523K**  
170 mm, 6 3/4"



**R7040D**  
**DU543K**  
135 mm, 5 1/4"




FREER  
R7050D

160 mm, 6 1/4"




FREER  
R7060D

160 mm, 6 1/4"




R7070D

145 mm, 5 3/4"




MEAD  
R7080D

165 mm, 6 1/4"




MOLT  
R7090D

180 mm, 7"




ROBERTS  
R7100D

185 mm, 7 1/4"




SELDIN  
R7110D

195 mm, 7 3/4"



SELDIN  
R7120D

195 mm, 7 3/4"




HOWARTH  
R7130D

215 mm, 8 1/2"



**LAMBOTTE  
R7140D**  
4 mm  
6 mm  
8 mm  
10 mm  
12 mm  
15 mm

Osteotome  
Osteotomes  
125 mm, 5"




**R7150D**  
4 mm  
6 mm  
8 mm  
10 mm  
12 mm

Osteotome  
Osteotomes  
140 mm, 5 1/2"




**PARTSCH  
R7160D**  
3 mm  
4 mm  
5 mm  
6 mm

Flachmeißel  
Chisels  
140 mm, 5 1/2"




**PARTSCH  
R7170D**  
3 mm  
4 mm  
5 mm  
6 mm  
7 mm

Chisels  
170 mm, 6 3/4"



**PARTSCH  
R7180D**  
3 mm  
4 mm  
5 mm  
6 mm

Gouges  
140 mm, 5 1/2"



**PARTSCH  
R7190D**  
3 mm  
4 mm  
5 mm  
6 mm  
7 mm

Gouges  
170 mm, 6 3/4"




**KHOURY**  
**R7200D** 4 mm  
6 mm

Chisel, curved on flat  
delicate model for pre-implantological operations  
125 mm, 5"




**KHOURY**  
**R7210D** 4 mm  
6 mm

Chisel, bajonet-shaped  
delicate model for pre-implantological operations  
especially in the posterior part.  
125 mm, 5"




**FREER**  
**R7220D** 4 mm

Chisel  
160 mm, 6 1/4"



**FREER**  
**R7230D** 4 mm

Chisel, curved  
160 mm, 6 1/4"



**LUCAS**  
**R7240D** 4 mm  
6 mm

Chisel  
160 mm, 6 1/4"



**LUCAS**  
**R7250D** 3 mm  
5 mm

Gouges  
160 mm, 6 1/4"

○ R7300D

○ R7290D

○ R7280D

$\frac{1}{16}$  ○ R7270D

R7260D

Tampers for bone grafts  
can also be used as haemostatic chisels  
160 mm, 6  $\frac{1}{2}$ "



R7310D

Tampers for bone grafts  
can also be used as haemostatic chisels  
120 mm, 4  $\frac{3}{4}$ "

□ R7330D

$\frac{1}{16}$  □ R7320D



R7340D

11 mm

Pterygomaxillary osteotome, curved, for anterior displacement of maxilla

235 mm, 9  $\frac{1}{4}$ "



R7350D

8 mm  
10 mm

Osteotomes, curved, for sagittal splitting procedure of the lower jaw  
250 mm, 10"



COTTLE  
R7360D

Osteotome with depth markings  
180 mm, 7"



COTTLE  
R7370D

Chisel-Osteotomes with depth markings  
180 mm, 7"



COTTLE  
R7380D

Osteotome with depth markings  
180 mm, 7"

Weights refer to the heads of the mallets




**R7390D**  
165 mm, 6 1/2"




**R7400D**  
190 mm, 7 1/2"



**R7410D**  
with interchangeable  
plastic disks  
185 mm, 7 1/4"  
Pair of spare disks



**MEAD**  
**R7420D**  
with interchangeable  
plastic disks  
190 mm, 7 1/2"



**COTTLE**  
**R7430D**  
1 side flat  
1 side convex  
185 mm, 7 1/4"




**HAJEK**  
**R7440D**  
220 mm, 8 3/4"

1/2 Pair of spare disks

1/2

1/2





**R7450D**  
Mallet,  
with interchangeable plastic disks  
185 mm, 7 1/4"

1/2 Pair of spare disks



**R7460D**  
Reposition forceps for adults  
for the left side,  
width of jaw 8,5 x 13 mm  
240 mm, 9 1/4"

1/2



**R7470D**  
Reposition forceps for adults  
for the right side,  
width of jaw 8,5 x 13 mm  
240 mm, 9 1/4"

1/2