


R1008D

Explorer,
delicate, flexible
155 mm, 6"




R1009D

Tooth tweezer,
serrated
160 mm, 6 1/4"



R1010D

Scaler for removing adhe-
sive and cement residues
175 mm, 7"



R1011D


Scaler/Pusher
175 mm, 7"

SCHURE



R1012D

straight
140 mm, 5 1/2"




R1013D

curved
140 mm, 5 1/2"



R1014D

angulated
120 mm, 5"



ADERER
R1015D

The snap device at the band seaters permits individual positioning of the bands



R1016D



R1017D



R1018D



R1019D

attachment DP793R
185 mm, 7 1/4"

2.5 x 1.5 mm

attachment DP794R
185 mm, 7 1/4"



R1020D

Separating forceps for inserting elastic separators
175 mm, 7"



R1021D

Ligature cutter with blades closing exactly for hard wire up to 4.0 mm diam.
115 mm, 4 1/2"



R1022D
Arch cutter


for hard wire up to 0.7 mm diam.
115 mm, 4 1/2"




ADERER

R1023D


Anterior band removing pliers for anteriors
130 mm, 5 1/8"



R1024D



TWEED
R1025D



JARABAK
R1026D

Wire bending pliers for shaping and closing loops with 3 holding grooves for wires up to 0.7 mm diam.
140 mm, 5 1/2"



ANGLE
R1027D


Bird-Beak-Pliers
Wire bending pliers for wire up to 0.7 mm diam.
125 mm, 5"



ADAMS
R1028D



ANGLE / TWEED
R1029D



HOW
R1030D

Wire fixation forceps straight, limb projecting on one side, universal forceps for multiband technique
145 mm, 5³/₄"



HOW
R1031D


Wire fixation forceps curved, limb projecting on one side, universal forceps for multiband technique
140 mm, 5¹/₂"

Universal forceps
serrated jaws, slim and curved for ideal positioning and distal bending of arches




WEINGART
R1032D

140 mm, 5¹/₂"



KORKHAUS
R1033D

Wire and ligature forceps
165 mm, 6 1/2"



R1034D

Ligature forceps, delicate, serrated jaws
for applying rubber ligatures
140 mm, 4 3/4"



MATHIEU
R1035D

Ligature forceps
140 mm, 5 1/2"



HALSTED-MOSQUITO
R1036D

Ligature forceps
125 mm, 5"




R1037D

Wire bending pliers,
serrated jaws
120 mm, 4 3/4"




MARBURG MODELL



R1039D

smooth jaws
140 mm, 5 1/2"




R1040D

jaws with
tungsten carbide inserts
140 mm, 5 1/2"


R1038D

Bird-Beak-Pliers
Wire bending pliers,
serrated jaws
140 mm, 5 1/2"



R1041D


Wire bending pliers
140 mm, 5 1/2"



GOSLEE

R1042D

Wire bending pliers, one smooth jaw,
one serrated jaw
140 mm, 5 1/2"



PEESO

R1043D

Wire bending pliers
135 mm, 5 1/4"



R1044D

Round nose pliers fine pattern
140 mm, 5 1/2"



LANGENBECK

R1045D


Round nose pliers
140 mm, 5 1/2"

Bending Pliers, Instruments for the Removable Technique



R1046D

Wire bending pliers
serrated jaws



R1047D
140 mm, 5 1/2"

R1048D
145 mm, 5 3/4"

R1049D
140 mm, 5 1/2"




ANDRESEN

1,0 mm Ø / 1.0 mm diam.

R1050D
Wire bending pliers
for bending loops
130 mm, 5 1/8"


Bending Pliers, Instruments for the Removable Technique



SCHWARZ

R1051D

Wire bending pliers
for bending loops
135 mm, 5 1/4"



FISCHER

R1052D

Wire bending pliers
for bending loops
130 mm, 5 1/8"



ADERER

R1053D

Wire bending pliers
for wire up to 0.7 mm diam
120 mm, 4 3/4"



ADERER

R1054D


Wire bending pliers, heavy pattern
for wire up to 0.9 mm diam
125 mm, 5"



ADERER

R1055D

Clasp forming and bar bending pliers
for oval wire up to 1.5 x 0.75 mm
145 mm, 5³/₄"



R1056D

Triangular bending clamp pliers
120 mm, 4³/₄"



YOUNG

R1057D

Clasp forming and bar bending pliers for hard wire up to 0.7 mm diam.
140 mm, 5¹/₂"



R1058D

Clasp forming pliers for springy wire up to 0.7 mm diam.
140 mm, 5¹/₂"

Bending Pliers, Instruments for the Removable Technique



R1059D

Pin- and wire cutters
for soft wire up to diam 1.5 mm, for hard wire up to diam 1.0 mm
150 mm, 6"




R1060D

Pin- and wire cutters
for soft wire up to diam 1.5 mm, for hard wire up to diam 1.2 mm
160 mm, 6 1/4"



R1070D

Arched pliers, for bending thicker wires (e.g. face bows)
200 mm, 8"



R1071D

Parallel flat pliers, for bending thicker wires (e.g. face bows)
185 mm, 7 1/4"

Wire Cutters, Instruments for the Removable Technique, Orthodontics Pliers



R1072D

Wire cutter
for soft wire up to diam 1.5 mm,
for hard wire up to diam 1.0 mm
145 mm, 5 $\frac{3}{4}$ "



R1073D

Wire cutter
for soft wire up to diam. 1.5 mm,
for hard wire up to diam. 1.0 mm
140 mm, 5 $\frac{1}{2}$ "



ANGLE
R1074D

Wire cutter
for soft wire up to 0.5 mm,
for hard wire up to 0.3 mm
130 mm, 5 $\frac{1}{8}$ "




ANGLE
R1075D


Wire cutter
for soft wire up to 0.6 mm,
for hard wire up to 0.4 mm
130 mm, 5 $\frac{1}{8}$ "



R1076D
120 mm, 4³/₄"



R1077D
Arrow clasp forming pliers



R1078D
Arrow clasp bending pliers
for round, hard wire up to 0.7 mm



ADAMS / NANCE

R1079D
Wire bending pliers
for hard wire up to 0.7 mm diam.
135 mm, 5¹/₄"



WALDSACHS


R1080D
Wire bending pliers
for wire up to 1.0 mm diam.
150 mm, 6"



PEESO

R1081D

Wire bending pliers
135 mm, 5 1/4"



TWEED


R1082D

Wire bending pliers for wire up to 0.9 mm diam.
135 mm, 5 1/4"



R1083D

Wire bending pliers
125 mm, 5"



R1084D


Wire bending pliers
135 mm, 5 1/4"



MC KELLOPS

R1085D

Wire bending pliers
160 mm, 6 1/4"



GORDON

R1086D

Contouring pliers
135 mm, 5 1/4"



REYNOLDS

R1087D

Contouring pliers
130 mm, 5"




JOHNSON

R1088D


Contouring pliers
135 mm, 5 1/4"

Wire and Ligature Forceps, Wire Scissors, Wire Cutters




R1089D

Wire scissors, 1 blade serrated
for soft wire up to diam. 1.2 mm,
for hard wire up to diam. 0.6 mm
120 mm, 4³/₄"



R1090D

Wire scissors, 1 blade serrated
for soft wire up to diam. 1.2 mm,
for hard wire up to diam. 0.8 mm
120 mm, 4³/₄"



TC
jaws with
tungsten carbide inserts



REILL

R1091D

Wire cutter TC
for soft wire up to diam. 2.0 mm
for hard wire up to diam. 1.5 mm
175 mm, 7"



R1092D

Wire cutter TC
for soft wire up to diam. 3.0 mm
for hard wire up to diam. 2.2 mm
235 mm, 9¹/₄"